

Correlation of fluorochromes and filter cubes

| Fluochrome | Filter cube |
|--|---------------|
| – Acid fuchsin | N 2.1, M 2 |
| – Acridine blue | A |
| – Acridine yellow | I 3, H 3 |
| – Acridine orange | I 3, H 3 |
| – Acridine red | N 2.1, N 3 |
| – Acriflavin | E 4, H 3 |
| – Acriflavin-Feulgen-SITS (AFS) | D |
| – Alizarin complexon | N 2.1 |
| – Alizarin red | N 2.1 |
| – Allophycocyanin (APC) | Y 3, Y 5 |
| – AMCA (Aminocoumarin) | A |
| – AMCA/FITC/Texas Red | B/G/R |
| – Aminoactinomycin D-AAD | N 2.1, N 3 |
| – Aniline blue | A |
| – ACMA | E 4 |
| – Astrazone Brilliant Red 4G | N 2.1 |
| – Astrazone Red 6B | N 2.1 |
| – Astrazone Yellow 7 CLL | H 3 |
| – Astrazone Orange R | I 3, L 5 |
| – Atabrine | E 4, H 3 |
| – Auramine | I 3, H 3 |
| – Aurophosphine, Aurophosphine G | I 3, H 3 |
| – BCECF | L 5 |
| – Berberine sulphate | H 3 |
| – Benzoxanthen Yellow | D |
| – BisAminophenyl Oxidiazol (BAO) | A |
| – Bisbenzimidazole (Hoechst) | A, D |
| – Blancophor BA | A, D, H 3 |
| – Blancophor SV | A |
| – BODIPY FL | L 5, K 3, I 3 |
| – Brilliant Sulphaflavine FF | D, H 3 |
| – Bromobimane (Thiolyte) | D |
| – Calcein | I 3 |
| – Calcein blue | A |
| – Calcium Crimson | Y 3 |
| – Calcium Green | K 3, I 3, L 5 |
| – Calcium Orange | M 2, N 2.1 |
| – Calcofluor White | H 3, D |
| – Calcofluor White standard solution | A |
| – Carboxyfluorescein diacetate C-FDA | I 3, L 5 |
| – Cascade Blue | A, D |
| – Catecholamines (adrenalin, noradrenalin, dopa, dopamine) | D |
| – Chromomycin A (mithramycin, olivomycin) | E 4 |
| – Coriphosphine O | I 3, H 3 |
| – Coumarin-phalloidin | D |
| – Cy 3 | Y 3 |
| – Cy 5 | Y 5 |
| – Cy 7 | Y 7 |

Correlation of fluorochromes and filter cubes

| Fluochrome | Filter cube |
|---|--------------------|
| – DANS (diamino-naphtyl sulphonic acid) | A |
| – DAPI | A, D |
| – DAPI (selective) | A 4 |
| – DAPI/FITC/Texas Red (simultaneous) | B/G/R |
| – Dansyl chloride | A |
| – DIPI | A |
| – DiI | Y 3 |
| – DiO | I 3, K 3 |
| – Diphenyl brilliant flavine 7 GFF | H 3 |
| – Dopamine | A |
| – DPH (diphenyl hexatriene) | A |
| – Eosin | B N 2.1 |
| – Ethidium bromide | N 2.1 |
| – Euchrysin | H 3, D |
| – Evans Blue | N 2.1 |
| – Fast Blue | A |
| – Fast Green FC G | N 2.1, M 2 |
| – Feulgen | N 2.1, TX 2 |
| – FDA (fluorescein diacetate) | I 3, H 3, K 3, L 5 |
| – FIF (formaldehyde induced fluorescence) | D, A |
| – FITC (fluorescein isothiocyanate) | I 3, H 3, K 3, L 5 |
| – FITC/ethidium bromide | I 3, L 5, N 2.1 |
| – FITC/phycoerithrin (PE) (simultaneous) | G/R |
| – FITC/Texas Red (simultaneous) | G/R |
| – FITC/TRITC (simultaneous) | FI/RH |
| – FITC (selective) | L 5 |
| – Texas Red (selective) | TX 2 |
| – FITC/TRITC | L 5, N 3 |
| – TRITC (selective) | N 3 |
| – Fluo 3 | I 3, L 5 |
| – Fluoro Gold | A |
| – Floram (fluorescamine) | A |
| – Genacryl Brilliant Red B | N 2.1 |
| – Genacryl Brilliant Yellow | E 4 |
| – Generic Blue | D |
| – GFP (Green Fluorescent Protein) | GFP |
| – Granular Blue | A |
| – Haematoporphyrin | N 2.1 |
| – Hoechst dye no. 33258 | A, D, A 4 |
| no. 33342 | A, D, A 4 |
| – Hydroxy-4-methylcoumarin | A |
| – Lissamine-rhodamine B (RB 200) | N 2.1, M 2 |
| – Lucifer Yellow | E 4 |
| – Magdala Red | N 2.1 |
| – Maleimide | A |
| – Mepacrin | D |
| – Merocyanin 540 | N 2.1 |
| – Mithramycin | E 4 |

Correlation of fluorochromes and filter cubes

| Fluochrome | Filter cube |
|--|-----------------|
| – MPS (methyl Green Pyronine stilbene) | A |
| – Nile Red | I 3, L 5, N .21 |
| – Nuclear Fast Red | N 2.1, M 2, N 3 |
| – Nuclear Yellow | A |
| – Olivomycin | E 4 |
| – Oregon Green (488, 500, 514) | L 5 |
| – Oxytetracycline | D |
| – Pararosaniline (Feulgen) | N 2.1, TX 2 |
| – Phosphine 3 R | I 3, H 3 |
| – Phycoerythrin (PE) | N 2.1, N 3 |
| – Primulin O | D |
| – Procion Yellow | D, E 4, H 3 |
| – Propidium iodide | N 2.1 |
| – Pyronine B | N 2.1, M 2 |
| – Quinacrine mustard (QM) | E 4 |
| – Resorufin | N 2.1, Y 3 |
| – Reverine | D |
| – Rhodamine B | N 2.1 |
| – Rhodamine 123 | I 3, L 5 |
| – Serotin | A, D |
| – SITS (stilbene isothiosulphonic acid) | A |
| – SITS acriflavine Feulgen | D |
| – Spectrum Orange | M 2, N 2.1 |
| – Sulphaflavine | A |
| – Tetracyclines: oxytetracycline, tetracycline, reverine (pyrrolidinomethyltetracycline), chlortetracycline, dimethylchlortetracycline | D |
| – Texas Red | TX 2 |
| – Thiazin red R | N 2.1, M 2 |
| – Thioflavine S | H 3, D |
| – Thioflavine TCN | A |
| – Thiolyte (bromobimane) | D, A |
| – TRITC (tetramethyl rhodamin isothiocyanate) | N 2.1, N 3 |
| – TRITC (selective) | N 3 |
| – True Blue | A |
| – Uranine B | H 3 |
| – Uvitex 2 B | A, D |
| – XRITC | N 2.1, N 3 |
| – Xylene orange | N 2.1, M 2 |