





MoticEasyScan Pro 24

MoticEasyScanner

Other FDA Scanner







MoticEasyScan Pro 24

New 12 MP camera combined with Motic 40x optics delivers imagery up to 0.13 um/pixel

MoticEasyScanner


Other FDA Scanner






MoticEasyScan Pro 24

More heavily saturated colors make it easier to see red stained cells in the sea of pink

MoticEasyScanner



Other FDA Scanner



MoticEasyScan Pro 24

Motic imagery is sharper and appears to have greater depth of field



MoticEasyScan Pro 24

- High-volume upgradability
- Enhance ROI detection Accuracy
- High-volume frozen section live

Want to send your slides to us and see the slide image quality yourself?

contact rdm@meyerinst.com or call 281-579-0342